

Biomechanical Comparison of Soft Tissue Fixation Strength of Coronet™ Soft Tissue Fixation System and Suture

CoNextions Research & Development

Objective

The purpose of this study was to compare the ultimate pull out strength of Coronet Soft Tissue Fixation System with that of a traditional suture repair in cadaveric tendons of the foot using a laterally applied load.

Background

The current standard of care uses knotted suture to achieve soft tissue fixation during tenodesis procedures. Accordingly, the most common source of failure for these types of repair is suture pulling through tendon (1). The Coronet Soft Tissue Fixation System was designed to address this limitation inherent in many tenodesis devices by using an electropolished stainless steel washer (8.7 mm outer diameter, 0.5 mm thickness) with six tines spaced equidistance around the circumference of the washer to fixate to the underlying soft tissue at the repair site. The soft tissue washer is secured to a 3.5 mm diameter PEEK bone anchor via a continuous loop of #2 suture (Figure 1). The device is provided pre-assembled to an inserter to allow for simultaneous placement of the bone anchor and soft tissue washer during implantation.




Figure 1: Coronet Implant

Sample Preparation

A total of 10 tendons were harvested from the feet of 4 adult cadaveric specimens. Peroneus Brevis (PB), Peroneus Longus (PL), Posterior Tibial (PT), and Anterior Tibialis (AT) tendons

Coronet



Suture




Figure 2: Tensile Test Set-up

were selected due to the propensity for these tendons to be reattached to bone in surgical settings. Harvested tendons were sectioned along their transverse axis creating proximal and distal tendon sections. To reduce the impact of variable tissue quality on the testing, one section of each tendon sample was used to test the Coronet, with the other section of the same tendon sample used to test suture.

To reduce the impact of variable bone quality and ensure the failure of the soft tissue connection during testing, synthetic test blocks were used to represent bone for the study (40 PCF density Sawbones, Pacific Laboratories, Inc).

Coronet Samples: Coronet samples were prepared by attaching the tendon sections to the test blocks following the device's Instructions for Use (IFU Ref. PK01024).

Suture Samples: A single stitch pattern commonly used in clinical practice was placed on the tendon using #2 suture (FiberWire® Suture, Arthrex, Inc) (2). Suture from the tendon sections was threaded down a Ø2.8mm through hole in the test block, wrapped around a steel dowel pin, and threaded back up through the hole in the test block where it was secured with 3 square knots creating a simple interrupted stitch.

Tensile Testing

Pull tests were completed using an Instron® test frame (Instron 3342 Series Universal Testing Device, Instron Corporation) Figure 2 shows the test setup for the Coronet and Suture samples. The tensile load was applied laterally to the simulated soft-tissue bone fixation in order to be representative of the loading conditions seen clinically following various soft tissue to bone repairs performed in the foot and ankle (2). The pull rate of the applied load was set at 12.5mm/s.

Results

The failure mode for all test samples was pull-through at the soft tissue. All Coronet specimens were stronger than the corresponding suture repairs in the same cadaveric tissue. The average load at failure for the Coronet specimens was 130.6N, the average repair strength for the suture repairs was 65.5N (Table 1).

Force	Coronet (n=10)	Suture (n=10)
Max (N)	182.0	75.4
Min (N)	79.4	48.4
Avg (N)	130.6	65.5

Table 1: Summary of Pull Test Results

Conclusion

Coronet Soft Tissue Fixation System holds tissue at an average maximum applied load almost 2 times greater than a standard of care suture technique under lateral loading conditions in cadaveric tendon samples.

References

1. Cummins, Craig A. et al. Mode of failure for rotator cuff repair with suture anchors identified at revision surgery. *Journal of Shoulder and Elbow Surgery* 2003;12(2):128-133.
2. Dawson, David M. et al. Modified Kidner Procedure Utilizing a Mitek Bone Anchor. *JFAS* 1998;37(2):115-121.

This document is intended solely for the use of healthcare professionals. Benchtop testing may not be indicative of clinical performance. The information presented is intended to demonstrate a CoNextions Inc. product. CoNextions does not practice medicine and recommends that surgeons be trained in the use of any particular product before using it in surgery. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. CoNextions does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. Individual results will vary based on health, weight, activity and other variables. Not all patients are candidates for this product and/or procedure. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your CoNextions representative if you have questions about the availability of CoNextions products in your area. For indications, contraindications, warnings, precautions, potential adverse effects and patient counselling information, see the package insert or contact your local representative. Visit www.conexionsmed.com for additional product information. CoNextions owns, uses, or has applied for the following trademarks: CoNextions, Coronet, Revolutionizing Tendon Repair. Sawbones is a trademark of Pacific Laboratories Inc. FiberWire is a trademark of Arthrex, Inc. Instron is a trademark of Illinois Tool Works Inc.

Manufactured by:
CoNextions Inc.
150 North Wright Brothers Drive,
Suite 560
Salt Lake City, Utah 84116, USA

Ordering Information/
Technical Information:
+1 (385) 351 1461
www.conexionsmed.com